







-E-Chair Instruction Page-

Icons on your e-chair Joystick:

 <p>“Joystick Hand controller” Tilting this forward, sideways or backwards controls the direction you can drive your e-chair</p>	 <p>“High / Low speed” This limits your maximum speed on whichever setting you are on</p>	 <p>“Speed Knob” Turning towards the turtle is slower Turning towards the rabbit is faster</p>	 <p>“Battery Indicator” Green bars or the right side is fully charged. Red or the left side is empty. When it gets to the middle, start to finish your journey.</p>
 <p>“Horn”</p>	 <p>“Power on” Turns your device on.</p>		



Speed Lower / Higher

The speed buttons can be either a knob or a series of orange bars, as the example image shows.


Getting on and off the machine- ALWAYS turn your machine off before getting on. This is to avoid accidentally touching the driving pad / Hand controller, and moving the machine.

Electric wheelchairs have two brakes at the back, found near each back tyre. If something goes wrong with the chair, take the E-Chair out of gear and push user back to the car. Gear Lever locations are found on your [specific rental device guides](https://www.mobilityscootersforhire.com.au/device-guides.php).

<https://www.mobilityscootersforhire.com.au/device-guides.php>

E-Chair Advice:

1. Lift the footrest/s up before getting on. They are NOT steps for the user to get onto the machine.
2. Drive slow by checking your **Speed Control Knob / Buttons** until you are familiar with using your e-chair.
3. When crossing roads/footpaths and going over bumps, always cross straight, so both your front wheels touch the road at the same time. This is to **avoid tipping your scooter**. **NEVER Go at an angle.**





CAUTION: This is ONLY for machines which CAN cross roads. Always adhere to the specs of your rental.

Driving Example B:
Bad! Your front left wheel hits the road first, tilting your scooter a bit. Then, your right front, and then your back left. Potentially can tip your machine sideways

-E-Chair Instruction Page-

4. When preparing to cross driveways and pedestrian footpaths, do not go too fast, and do NOT go to slow.
 - a. If you stop right before going over the curb, your machine won't have enough power to get you up.
 - b. If you go too fast when going over the curb, you might lose control of your tiller.
 - c. **Always go at a steady pace that can be controlled.**
5. If your e-Chair can take hills / inclines or you have a steep driveway to get over, **put the speed down when going upwards on these slopes.** The e-chair draws a lot of power from the battery when going up these slopes, and to avoid stressing / wrecking the battery and overworking the machine, go at a lower speed. Your machine **will** slow down when going upwards, but this is normal. NEVER go at full throttle up these hills / inclines.
6. You cannot drive the e-chair if the charger is in. Please ensure it is out before turning it on.
7. Charge your e-chair after each trip. Always be aware of what type of scooter you have. AGM / GEL Deep cycle batteries mean you can leave your scooter on charge continuously. If your machine has **lithium Ion Batteries**, as soon as the charger has turned green, take your scooter off charge immediately.
8. If your machine has lithium Ion batteries, it will be found in the description of each rental available.

 <p>Faldan</p>	<p>Faldan a Compact Light Weight Folding Electric Power Chair. 28.5 kg total weight, 23.5 kg without batteries. Very compact and very maneuverable allowing the driver to handle tightest of environments inside, or outside, on smooth flat surfaces. Faldan folds easily and quickly for transport or storage. Two Li-Ion batteries certified for airplane use, making the Faldan truly a great transportable. 2x 6Ah Li-Ion. 6 kph. 130 kg load. Please note you cannot put a rear bag on this one.</p>	<p>Make a booking</p> <p>26/11-11/02+</p>	<p>75.00</p>
 <p>Blueberry</p>	<p>Power Electric Wheelchairs™ Office chair on wheels! Front wheel drive, compact, comfortable office chair with seat belt, designed for inside / outside use, able to tackle slopes. 35 Ah batteries, 6 kph, 115 kg load</p>	<p>Make a booking</p>	<p>70.00</p>

Batteries	1x 25V 1.06 Ah 274 Wh Li-Ion (Std on all models)
Batteries	2 x 12 V x 18 Ahr

You can see how this machines clearly state Lithium ion batteries. Lithium ion Batteries can also be written as "Li-Ion Battery".

You can see how this machine just say 12Ah batteries, so they are not lithium.

If you are ever unsure, you can also click on the image or the name of the e-chair and have a look in the specifications at the bottom of the new page. **Alternatively, just ask us or simply take the machine off charge when its fully charged anyway.**

9. It is **VERY** important to check, if you are going on a cruise ship, to ensure the volt range is **110V-240V**, so your chargers can turn on.
10. **Do not leave your e-chair unattended.** It is your responsibility to pay for a replacement e-chair, should it get lost or stolen. Also, **ensure you maintain contactable with Mobility Scooters For Hire** for the duration of your hire, should we need to contact you, or visa versa.
11. **If you crash, and parts break off, please pick the pieces up.** We may be able to repair. Your Insurance covers you for 65% of damage repair, but NOT replacement of parts you left behind.
12. It is **recommended** you download the **E-Chair Guide** for your machine as well, so you have advice beyond the instructions on this sheet and the demonstration our Team gives you, should we be outside business hours. Go to the webpage and select your e-Chair:
<https://www.mobilityscootersforhire.com.au/device-guides.php>
13. The name of your e-chair is on the Joystick or side panels of your e-chair, on your battery box for some portable machines, and always on your invoice. When speaking to an office attendant about your e-chair, please refer to the e-chair name or the name of the person the invoice is booked under. Alternatively, you can provide the invoice number.
14. If you **require an extension**, please contact us 1 day prior to the end date. Your booking dates are on your invoice, so please keep track of your dates.
 - a. Please note we will usually charge your credit card for payment of the extension invoice, so please ensure you have funds available on the day of extension if you have asked for an extension.