

# -Knee Scoot Instruction Page-



## Knee Scoots Instructions:

1. **Do not leave your Knee Scoot unattended.** It is your responsibility to pay for a replacement, should it get lost or stolen.
2. **Knee Scoots are designed to be used ONLY on flat ground.**
  - a. Any change in ground level, whether it be a "lip" of only 3 or 4mm between the kitchen and lounge room, or going over a carpet or power cord on the floor of the house, stop, and carefully go over it with your machine.
  - b. Do not use it on footpaths outside, as it is not designed for outdoor use. The footpaths are uneven and some parts sink or go higher.
  - c. **LOAD RESTRICTIONS** – The load in the front basket must not exceed 1kg.
3. **Using it on slopes - Your Knee Scoot brakes are NOT designed to stop you crashing or to slow down your speed.** They are only for when you are getting on/off the machine, so your Knee Scoot does not move. As such, you should always go at a walking pace, so you can control your machine.
  - a. You may require a disability ramp to get into your doctor's clinic or hospital. In these cases, it is only a few metres, and all smooth flat ground, with a railing beside you. You can use your Knee Scoot here. Always go up or down carefully, at a walking pace.
4. **Always have both hands on the steering tiller,** so you have complete control over the knee scoot.
5. To stop your Knee scoot, you should use a combination of your handbrakes AND your foot touching the ground to slow yourself. You would only be at a walking pace, and no faster.
6. **Maintenance** - Sometimes you may have the knee scoot for a longer than expected timeframe. Ensure you tighten any loose nuts or bolts, or alternatively organise with us a time to bring it down for a quick service.
7. **Proper Shoe Wear Equipment** – The foot you use to kick off from the ground and helps you "scoot" along : you should wear runners or a form of shoe that has grip or traction on the ground. Do not wear shoes that will make it harder for you to stop, should the need arise.
8. **Your Knee Scoot is not designed for sitting** – the kneepad area is ONLY for you to place your knee. knee scooter is designed to carry ONE person only.
9. **To Fold your Knee Scoot:**

### 2. Pull the lever out



1. **Lift the Silver Clasp Up.**



3. **Pull to the left / right to release the pin, and fold tiller down.**

