

4-WHEEL SCOOTER

# PATHRIDER 130 XL

OUTSTANDING OUTDOOR  
PERFORMANCE WITH  
LUXURIOUS FEATURES

The Pathrider® 130 XL features a powerful drivetrain and large 33 cm (13") pneumatic wheels to handle rugged outdoor terrain. Luxury touches like a wraparound easy drive tiller, full suspension and seat post suspension for a smoother ride and an LED kerb light blend seamlessly giving the Pathrider 130 XL a sleek, sporty style.



33 cm (13")  
Pneumatic Tyres



Full Suspension



Easy Drive Tiller



Full lighting package

PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

[www.pridemobility.com.au](http://www.pridemobility.com.au)

# SPECIFICATIONS

## PATHRIDER® 130 XL

## 4-wheel

Weight capacity:	158.75 kg (350 lbs.)
Maximum speed <sup>1</sup> :	Up to 15 km/h (9.5 mph)
Ground clearance <sup>2</sup> :	9.69 cm (3.81")
Turning radius <sup>2</sup> :	170 cm (67")
Overall length <sup>2</sup> :	136 cm (53.54")
Overall width <sup>2</sup> :	63.65 cm (25")
Seat-to-floor height:	63 cm - 67 cm (25" - 26")
Front tyres:	10 cm x 33 cm (4" x 13") pneumatic
Rear tyres:	10 cm x 33 cm (4" x 13") pneumatic
Suspension:	Front and rear
Range per charge <sup>1,4</sup> :	Up to 39.5 km (24.58 miles)
Maximum allowable incline:	10°
Total weight without batteries <sup>2,5</sup> :	65.28 kg (143 lbs.)
Heaviest piece when disassembled:	43 kg (94.5 lbs.) (front section)
Standard seating:	Type: Pillow top deluxe high-back with sliders Weight: 21.77 kg (48 lbs.) Material: Black Vinyl Dimensions: 49.5 cm (19.5") width, 44.5 cm (17.5") depth
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	(2) 12-volt, deep cycle Size: 55 AH (standard) Weight <sup>6</sup> : 16.6 kg (36.5 lbs.) each
Battery charger:	5 amp, off-board charger
Colours:	



Mercury Silver



Candy Apple Red

## FEATURES

- Large 33 cm (13") pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, dual bulb lower LED headlight
- Seat post suspension for a smoother ride
- Easy drive tiller with wraparound handles
- Front basket standard
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror standard

## OPTIONS

- Rear basket
- Safety flag
- Walker holder

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+10%, -5%).  
<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).  
<sup>3</sup> AGM or gel cell type required.  
<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.  
<sup>5</sup> Includes standard seat.  
<sup>6</sup> Battery weight may vary based on manufacturer.

**NOTE: All specifications subject to change without notice.**  
 The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®



(Tel) 03 8770 9600 • (Fax) 03 9703 2720

20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

[www.pridemobility.com.au](http://www.pridemobility.com.au)

