

Invacare® Auriga Scooter



Yes, you can.®

Auriga Scooter





Yes, you can.®

Invacare Canada
www.invacare.ca

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324



3-Wheel Version: S8301-ARIN - Red
S8301-SBIN - Blue



4-Wheel Version: S8401-ARIN - Red
S8401-SBIN - Blue

The Invacare® Auriga is highly maneuverable, offering spacious legroom, a powerful motor and a maximum user weight of 330 lbs. The Auriga ensures a safe driving experience with controls that are adjustable and easy to read and buttons that are easy to press and respond immediately. With many adjustable options to meet your requirements, your driving position is completely adjustable for your safety, comfort and for easy transfer. The Auriga is available in a 3-wheel version and a 4-wheel version. Designed for comfort and safety, the Auriga allows you daily independence in style.

Features

- Revolving Seat & Swing-up Armrest - enables easy access getting in and out of the seat
- Reclining Backrest - adjustment for individual comfort
- Adjustability - optimal driving position can be obtained by adjusting seat height, horizontal seat position and the width and height of the armrests
- Tiller - a unique design of the tiller allows the user to easily adjust the tiller angle to the optimal driving position
- Charger Status Window - a green light indicates "charging"
- Easy Transport - a quick release mechanism easily dismantles the scooter to fit into a car, allowing the scooter to be transported
- Handle Bar - an ergonomic handle bar enables light and easy turning
- Throttle - can be operated by a light "pull action"
- Safety Belt - provides support while driving
- The Auriga 4-Wheel version also incorporates an adjustable spring suspension which ensures smooth, comfortable driving with a top speed of 10 km/h
- Fully programmable for individual client needs



TILLER

Independent of the Tiller angle the scooter headlight always shines directly ahead.



CONTROL PANEL

The spacious control panel includes:

- "Easy touch" buttons
- Illuminating indicators when button is activated
- Accurate easy to read battery indicator
- Handy "Speed dial"



SEAT POSITIONING

7" sliding rails allows adjustment of the seat position.



ADJUSTABLE ARMRESTS


The "flip-up" armrests are width, angle and height adjustable.

Specifications

	Auriga 3-Wheel	Auriga 4-Wheel
Width	24"	24"
Height	47"	47"
Length (w/o bumpers and anti-tippers)	47"	48"
Seat width	17-21 1/2"	17-21 1/2"
Seat depth	17 1/2"	17 1/2"
Seat height (measured from the chassis, manually adjustable)	17, 18, 19"	17, 18, 19"
Backrest height (w/o headrest)	17"	17"
Backrest height (with headrest)	25"	25"
Backrest angle (manually adjustable)	0°, 10°, 15°, 30°	0°, 10°, 15°, 30°
Armrest height	9-11"	9-11"
Weight, including batteries (lbs.)	194	208
Maximum weight capacity (lbs.)	330	330
Speed	10 km/h	10 km/h
Batteries	2 x 40 Ah	2 x 40 Ah
Tire pressure	40 psi	40 psi
Minimum turning radius	47"	59"
Safe climbing capability 265 lb load	12°	12°
Safe climbing capability up to 330 lb load	10°	10°
Maximum obstacle height	3"	3"

© 2006
Invacare Canada

Form No. 05-450C



"Invacare says 'Yes, you can'.
and I couldn't have said it
better myself."