

**Things to do. Places to go.
Trust Invacare to take you.**



SCOOTER SPECIFICATION CHART

**FREE 24 HOUR
ROADSIDE
ASSISTANCE
PROGRAM***



Invacare Australia's stylish mobility scooter range, is internationally recognised for its definitive blend of safety, comfort, solid trusted design and secure manoeuvrability. Whether you are looking for a small, low-cost mobility scooter for short drives and getting around the shops, or you need a larger, higher-powered scooter that's built to suit the varied Australian terrain, the Invacare range of scooters has it all.




Every Invacare scooter is backed up by an efficient service network and FREE 24-hour roadside assistance for 14 months*! So when there are things to do and places to go, you can trust Invacare for all the worry-free mobility you want.






Yes, you can®

For further information just ask your Invacare Dealer, or visit our website www.invacare.com.au

* Conditions apply.

	Heavy Duty		Large
			
	Meteor	Meteor Hill Climber	Comet
Model	4 wheel	4 wheel high torque motor	4 wheel
Colours	Sapphire Blue, Ruby Red	Sapphire blue, Ruby Red	Moonstone Silver, Onyx Blue
Overall Length	146cm	146cm	147.5cm
Overall Width	66cm	66cm	66cm
Wheels: Front & Rear	31cm/ 12" Pneumatic	31cm/ 12" Pneumatic	33cm/ 13" Pneumatic
Weight without batteries	85kg	85kg	84kg
Max. Speed	15 km/h	10 km/h	15 km/h
Weight Capacity	200kg	200kg	160kg
Ground Clearance	11cm	11cm	13.5cm
Grade Climbable	10°	15°	10°
Kerb Climbing	11cm	11cm	10cm
Turning Radius	161cm	161cm	163cm
Suspension	Front & Rear, Spring & Adjustable	Front & Rear, Spring & Adjustable	Front & Rear, Spring & Adjustable
Seat Type	High back with Headrest	High back with Headrest	Extra soft, ergonomic, recline lever, high back with Headrest
Seat Width	44cm	44cm	51cm
Armrest Width	44-55cm	44-55cm	51-59cm
Continuous & Peak Motor Size	Cont. -400W, Peak - 1750W	Cont. -400W, Peak - 1750W	Cont. -250W, Peak - 1200W
Standard Battery Size	(2) 12V 73AH	(2) 12V 73AH	(2) 12V 75AH
Standard Travel Range	Up to 35km	Up to 35km	Up to 55km
Standard Battery Charger	DC 24V / 8 Amp	DC 24V / 8 Amp	DC 24V / 8 Amp
Electronics	Rhino 160 Amp Controller	Rhino 160 Amp Controller	Rhino 110 Amp Controller
Limited Warranty	Scooters 2 Years, 13 months MK Gel Batteries, Charger 1 year	Scooters 2 Years, 13 months MK Gel Batteries, Charger 1 year	Scooters 2 Years, Batteries & Charger 1 year

Medium		Small
		
Auriga	Pegasus	Leo
3 or 4 wheel	4 wheel	4 wheel
Sapphire blue, Ruby Red	Moonstone Silver, Onyx Blue	Moonstone Silver, Sapphire Blue, Ruby Red
3 Wheel: 120cm; 4 Wheel: 122cm	129cm	122cm
61cm	62.5cm	59cm
26cm/ 10" Pneumatic	29cm/ 11" Pneumatic	26cm/ 10" Pneumatic
3 Wheel: 66kg; 4 Wheel: 72kg	78kg	59kg
3 Wheel: 8km/h; 4 Wheel: 10km/h	10km/h	8km/h
150kg	136kg	136kg
10cm	9cm	7cm
3 Wheel: 10°; 4 Wheel: 12°	12°	10°
8cm	8cm	6cm
3 Wheel: 120cm; 4 Wheel: 150cm	140cm	131cm
Rear Spring & Adjustable	Front spring; Rear spring & Adjustable	N/A
Mid-back with Headrest	Extra soft, ergonomic, recline lever, high back with Headrest	Ergonomic, mid-back seat
44cm	51cm	45cm
44-55cm	51-59cm	45-55cm
Cont. -230W, Peak - 715W	Cont. -240W, Peak - 600W	Cont. -240W, Peak - 450W
(2) 12V 40AH	(2) 12V 50AH	(2) 12V 36AH
3 Wheel: Up to 40km; 4 Wheel: Up to 32km	Up to 42km	Up to 36km
DC 24V / 8 Amp	DC 24V / 8 Amp	DC 24V / 8 Amp
3 Wheel: Rhino 70 Amp Controller; 4 Wheel: Rhino 110 Amp Controller	Rhino 110 Amp Controller	Dynamic DR-90 90Amp Controller
Scooters 2 Years, 13 months MK Gel Batteries, Charger 1 year	Scooters 2 Years, Batteries & Charger 1 year	Scooters 18 months, Batteries & Charger 1 year



Yes, you can[®]

Invacare New Zealand
4 Westfield Place (PO Box 62-124) Mt Wellington Auckland, New Zealand
Freephone: 0508 468 222 Freefax: 0508 807 788
Email: sales@invacare.co.nz Web: www.invacare.co.nz

Invacare Australia Pty Ltd
1 Lenton Place (PO Box 5002) North Rocks NSW 2151, Australia
Freephone: 1800 460 460 Fax: 02 8839 5353
Email: sales@invacare.com.au Web: www.invacare.com.au

Invacare Asia Ltd
1 Lenton Place (PO Box 5002) North Rocks NSW 2151, Australia
Phone: +61 2 8839 5330 Fax: +61 2 8839 5331
Email: asiasales@invacare.com Web: www.invacare.com

All information is subject to change without prior notice. Photographs are for illustration purposes only and may vary from product description.
Invacare[®] is a registered trademark and Yes, you can.[™] is a trademark of Invacare[®] Corporation. 2009 Invacare Asia Pacific